

Imperial College
London

Nucleic Acid Aptamers for Sensing and Imaging


NOVEL DIAGNOSTICS FOR INFECTIOUS DISEASES
25-26 March, 2024; Royal Society, London

Tony Cass. Department of Chemistry

Aptamers-Strengths as Sensing Molecules

- Well-defined at the molecular level
 - Available in high quantity and quality via chemical synthesis
 - Low Cost (<1% that of antibodies per mole at a research scale)
 - Precision chemical modification either during, or post-synthesis
 - High stability (with suitable modification)
 - Straightforward Rules for Intramolecular Interactions
-

Aptamer Selection via SELEX




Other Flavours of SELEX:

1. Capture SELEX
2. Mag Bead SELEX
3. Microfluidic SELEX
4. CE-SELEX
5. Proximity Ligation SELEX

Labelling


Redox active groups
Fluorophores
Biotin
Nanoparticles

Truncation



5' Truncations

Fusions (Chimaeras)




Streptavidin~Theophylline
Binding RNA Aptamers


+ Theophylline

Aptamers in Virus Detection

3, or 5' Thiol Terminated Aptamer



Assay Design




Virus Specific Aptamer Random DNA Sequence



After George Whitesides

DRELFA: Dual Recognition Element Lateral Flow Assays



Standard Lateral Flow Format

1. Biotinylated aptamer for specificity
2. Streptavidin Capture Line
3. AuNP labelled antibody for detection
4. *anti* mouse antibody for QA
5. Avoids immobilised aptamer 'lines'

DRELFA Performance




(a)

Conventional 2-antibody






(b)

DRELFA




Aptamer *anti* H3N2/Panama
Antibody *anti* H3

Conformational Switching. Electrochemical Aptasensors




Imperial College London Aptamer Stability

G-Quadruplexes



G-quadruplex that binds
Retinoblastoma protein



Acceptor Photobleaching FRET

→ Intracellular
Stability

Aptamer ‘Specificity’ A Cautionary Tale

Nucleic Acid Aptamers are often claimed to be ‘Highly Specific’

Some Counter Examples

“Cocaine Binding Aptamer” **Binds quinine and cholesterol**

Nucleolin Binding Aptamer **A ‘control’ aptamer with a randomized sequence (AS1411) has a higher affinity**

Arsenite/arsenate Binding Aptamer **Binds neither, flawed assay**

Acknowledgements

Group: Thao Le, Orada Chumphukan, Kittichan Kanok, Shaolin Liang

Collaborators:

Andrew Ellington (UT Austin)

Munir Iqbal (Pirbright),

Waljit Dhillon, Liming Ying (ICL)

Julian Tanner (HKUST)
